

Product datasheet for **TR503679**

Ntmt1 Mouse shRNA Plasmid (Locus ID 66617)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Ntmt1 Mouse shRNA Plasmid (Locus ID 66617)
Locus ID:	66617
Synonyms:	2610205E22Rik; AL033331; AL033332; Mettl11a; NTM1A
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Ntmt1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 66617). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC027220 , NM_170592 , NM_001356433 , NM_001356434 , NM_001356435 , NM_001356436 , NM_001356437 , NR_151467 , NR_151468 , NM_170592.1 , NM_170592.2
UniProt ID:	Q8R2U4
Summary:	Distributive alpha-N-methyltransferase that methylates the N-terminus of target proteins containing the N-terminal motif [Ala/Gly/Pro/Ser]-Pro-Lys when the initiator Met is cleaved. Specifically catalyzes mono-, di- or tri-methylation of the exposed alpha-amino group of the Ala, Gly or Ser residue in the [Ala/Gly/Ser]-Pro-Lys motif and mono- or di-methylation of Pro in the Pro-Pro-Lys motif (PubMed:20668449). Some of the substrates may be primed by METTL11B-mediated monomethylation. Catalyzes the trimethylation of the N-terminal Gly in CENPA (after removal of Met-1) (By similarity). Responsible for the N-terminal methylation of KLHL31, MYL2, MYL3, RB1, RCC1, RPL23A and SET. Required during mitosis for normal bipolar spindle formation and chromosome segregation via its action on RCC1 (PubMed:20668449). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).