

Product datasheet for **TR503925**

Nudt7 Mouse shRNA Plasmid (Locus ID 67528)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Nudt7 Mouse shRNA Plasmid (Locus ID 67528)
Locus ID:	67528
Synonyms:	1300007B24Rik; 2210404C19Rik
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Nudt7 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 67528). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<u>BC069843</u> , <u>NM_001290180</u> , <u>NM_001290181</u> , <u>NM_001290182</u> , <u>NM_024437</u> , <u>NM_024446</u> , <u>NM_024437.1</u> , <u>NM_024437.2</u> , <u>NM_024437.3</u> , <u>NM_024437.4</u> , <u>NM_024446.1</u> , <u>NM_024446.2</u> , <u>NM_024446.3</u> , <u>NM_024446.4</u> , <u>NM_024446.5</u> , <u>NM_001290182.1</u> , <u>NM_001290180.1</u> , <u>NM_001290181.1</u> , <u>BC033046</u>
UniProt ID:	<u>Q99P30</u>
Summary:	Coenzyme A diphosphatase which mediates the cleavage of CoA, CoA esters and oxidized CoA with similar efficiencies, yielding 3',5'-ADP and the corresponding 4'-phosphopantetheine derivative as products. CoA into 3',5'-ADP and 4'-phosphopantetheine. Has no activity toward NDP-sugars, CDP-alcohols, (deoxy)nucleoside 5'-triphosphates, nucleoside 5'-di or monophosphates, diadenosine polyphosphates, NAD, NADH, NADP, NADPH or thymidine-5'-monophospho-p-nitrophenyl ester. May be required to eliminate oxidized CoA from peroxisomes, or regulate CoA and acyl-CoA levels in this organelle in response to metabolic demand. Does not play a role in U8 snoRNA decapping activity. Binds U8 snoRNA. [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).