

## Product datasheet for **TR504511**

### Tex19.2 Mouse shRNA Plasmid (Locus ID 70956)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Tex19.2 Mouse shRNA Plasmid (Locus ID 70956)
Locus ID:	70956
Synonyms:	4921530G04Rik; Tex19b
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Tex19.2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 70956). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_027622</a> , <a href="#">NM_027622.1</a> , <a href="#">NM_027622.2</a> , <a href="#">NM_027622.3</a> , <a href="#">BC141309</a>
UniProt ID:	<a href="#">Q9D5S1</a>
Summary:	May be required during spermatogenesis, probably by participating in the repression of retrotransposable elements and prevent their mobilization (Probable). With its paralog, Tex19.1, collaborates with the Piwi-interacting RNA (piRNA) pathway, which mediates the repression of transposable elements during meiosis by forming complexes composed of piRNAs and Piwi proteins. Interacts with Piwi proteins and directly binds piRNAs, a class of 24 to 30 nucleotide RNAs that are generated by a Dicer-independent mechanism and are primarily derived from transposons and other repeated sequence elements (PubMed:28254886).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).