

## Product datasheet for **TR505048**

### **Bcdin3d Mouse shRNA Plasmid (Locus ID 75284)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	Bcdin3d Mouse shRNA Plasmid (Locus ID 75284)
Locus ID:	75284
Synonyms:	4930556P03Rik; AV138748
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Bcdin3d - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 75284). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC016225</a> , <a href="#">NM_029236</a> , <a href="#">NM_029236.1</a> , <a href="#">NM_029236.2</a> , <a href="#">BC078455</a>
UniProt ID:	<a href="#">Q91YP1</a>
Summary:	O-methyltransferase that specifically monomethylates 5'-monophosphate of cytoplasmic histidyl tRNA, acting as a capping enzyme. Less efficiently, also methylates the 5' monophosphate of pre-miRNAs, acting as a negative regulator of miRNA processing. The 5' monophosphate of pre-miRNAs is recognized by DICER1 and is required for pre-miRNAs processing: methylation at this position reduces the processing of pre-miRNAs by DICER1. Able to mediate methylation of pre-miR-145, as well as other pre-miRNAs. There is some controversy about the methylation of pre-miR-145, since the dimethylation first described as the specific enzymatic activity cannot be reproduced by a more recent work which observes a monomehtylation of pre-miR-145 but two orders weaker than the methylation of cytosolic histidyl tRNA.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).