

Product datasheet for **TR505159**

Daam2 Mouse shRNA Plasmid (Locus ID 76441)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Daam2 Mouse shRNA Plasmid (Locus ID 76441)
Locus ID:	76441
Synonyms:	2310016D11Rik; AI843643; AW557870
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Daam2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 76441). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001008231 , NM_001008231.1 , NM_001008231.2 , BC037720 , BC050043
UniProt ID:	Q80U19
Summary:	Key regulator of the Wnt signaling pathway, which is required for various processes during development, such as dorsal patterning, determination of left/right symmetry or myelination in the central nervous system (PubMed:22227309, PubMed:24091014, PubMed:25754822). Acts downstream of Wnt ligands and upstream of beta-catenin (CTNNB1) (PubMed:22227309, PubMed:25754822). Required for canonical Wnt signaling pathway during patterning in the dorsal spinal cord by promoting the aggregation of Disheveled (Dvl) complexes, thereby clustering and formation of Wnt receptor signalosomes and potentiating Wnt activity (PubMed:22227309). During dorsal patterning of the spinal cord, inhibits oligodendrocytes differentiation via interaction with PIP5K1A (PubMed:25754822). Also regulates non-canonical Wnt signaling pathway (PubMed:24091014). Acts downstream of PITX2 in the developing gut and is required for left/right asymmetry within dorsal mesentery: affects mesenchymal condensation by lengthening cadherin-based junctions through WNT5A and non-canonical Wnt signaling, inducing polarized condensation in the left dorsal mesentery necessary to initiate gut rotation (PubMed:24091014). Together with DAAM1, required for myocardial maturation and sarcomere assembly (PubMed:26526197).[UniProtKB/Swiss-Prot Function]



[View online »](#)

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).