

Product datasheet for **TR505567**

Ssx2ip Mouse shRNA Plasmid (Locus ID 99167)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Ssx2ip Mouse shRNA Plasmid (Locus ID 99167)
Locus ID:	99167
Synonyms:	Adip; AU014939; AU042321
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Ssx2ip - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 99167). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC021749 , BC031527 , NM_001253768 , NM_001253769 , NM_001253770 , NM_138744 , NM_001355661 , NM_138744.1 , NM_138744.2 , NM_138744.3 , NM_001253768.1 , NM_001253769.1 , NM_001253770.1 , NM_001253770.2 , NM_001253768.2 , NM_001253769.2
UniProt ID:	Q8VC66
Summary:	Belongs to an adhesion system, which plays a role in the organization of homotypic, interneuronal and heterotypic cell-cell adherens junctions (AJs). May connect the nectin-afadin and E-cadherin-catenin system through alpha-actinin and may be involved in organization of the actin cytoskeleton at AJs through afadin and alpha-actinin (PubMed:12446711). Acts as a centrosome maturation factor, probably by maintaining the integrity of the pericentriolar material and proper microtubule nucleation at mitotic spindle poles. The function seems to implicate at least in part WRAP73; the SSX2IP:WRAP73 complex is proposed to act as regulator of spindle anchoring at the mitotic centrosome (By similarity). Involved in cell movement: localizes at the leading edge of moving cells in response to PDGF and is required for the formation of the leading edge and the promotion of cell movement, possibly via activation of Rac signaling (PubMed:22027834). Involved in ciliogenesis (By similarity). It is required for targeted recruitment of the BBSome, CEP290, RAB8, and SSTR3 to the cilia (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>