

Product datasheet for **TR506608**

Pcdh20 Mouse shRNA Plasmid (Locus ID 219257)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Pcdh20 Mouse shRNA Plasmid (Locus ID 219257)
Locus ID:	219257
Synonyms:	C630015B17Rik; PCDH1; Pcdh13
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Pcdh20 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 219257). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_178685 , NM_178685.1 , NM_178685.2 , NM_178685.3 , NM_178685.4 , NM_178685.5 , BC150665 , BC024927 , BC079605
UniProt ID:	Q8BIZ0
Summary:	This gene belongs to the protocadherin gene family, a subfamily of the cadherin superfamily. The encoded protein contains six extracellular cadherin domains, a transmembrane domain, and a cytoplasmic tail differing from those of the classical cadherins. The encoded protein may play a role in cell adhesion in the nervous system and has been shown to be specifically expressed in newly differentiated olfactory sensory neurons and their axons during development. In adult mice, the expression of this protein in the olfactory system is more restricted but shows a gender difference with higher expression in the male than in the female. [provided by RefSeq, Sep 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).