

Product datasheet for **TR506742**

Tiprl Mouse shRNA Plasmid (Locus ID 226591)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Tiprl Mouse shRNA Plasmid (Locus ID 226591)
Locus ID:	226591
Synonyms:	1810011K17Rik
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Tiprl - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 226591). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC002098 , BC058250 , NM_145513 , NM_145513.1 , NM_145513.2 , NM_145513.3 , NM_145513.4 , NM_001357410 , NM_001357411
UniProt ID:	Q8BH58
Summary:	May be a allosteric regulator of serine/threonine-protein phosphatase 2A (PP2A). Inhibits catalytic activity of the PP2A(D) core complex in vitro. The PP2A(C):TIPRL complex does not show phosphatase activity. Acts as negative regulator of serine/threonine-protein phosphatase 4 probably by inhibiting the formation of the active PPP4C:PPP4R2 complex; the function is proposed to implicate it in DNA damage response by promoting H2AFX phosphorylated on Ser-140 (gamma-H2AFX). May play a role in the regulation of ATM/ATR signaling pathway controlling DNA replication and repair (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).