

Product datasheet for **TR506772**

Dclre1c Mouse shRNA Plasmid (Locus ID 227525)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Dclre1c Mouse shRNA Plasmid (Locus ID 227525)
Locus ID:	227525
Synonyms:	9930121L06Rik; A; AI661365; Art; Snm1l
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Dclre1c - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 227525). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC108935 , NM_001110214 , NM_001302674 , NM_001302684 , NM_146114 , NM_175683 , NM_175683.1 , NM_175683.2 , NM_175683.3 , NM_175683.4 , NM_001110214.1 , NM_146114.1 , NM_146114.2 , NM_146114.3 , NM_001302674.1 , NM_001302684.1 , BC030672 , BC030680
UniProt ID:	Q8K4J0
Summary:	This gene encodes a member of the SNM1 family of nucleases and is involved in V(D)J recombination and DNA repair. This protein has single-strand-specific 5'-3' exonuclease activity; it also exhibits endonuclease activity on 5' and 3' overhangs and hairpins. The protein also functions in the regulation of the cell cycle in response to DNA damage. Homozygous knockout mice for this gene exhibit severe combined immunodeficiency with sensitivity to ionizing radiation. Mutations in this gene in humans can cause Athabaskan-type severe combined immunodeficiency (SCIDA) and Omenn syndrome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).