

## Product datasheet for **TR507132**

### Cplx3 Mouse shRNA Plasmid (Locus ID 235415)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cplx3 Mouse shRNA Plasmid (Locus ID 235415)
Locus ID:	235415
Synonyms:	BC016632; CpxIII; ERGIC-53L; ERGL; Lamn1I; Lman1I
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Cplx3 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 235415). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC016632</a> , <a href="#">BC024854</a> , <a href="#">BC033343</a> , <a href="#">NM_146223</a> , <a href="#">NM_146223.1</a> , <a href="#">NM_146223.2</a> , <a href="#">NM_146223.3</a> , <a href="#">BC017128</a>
UniProt ID:	<a href="#">Q8R1B5</a>
Summary:	Complexin that regulates SNARE protein complex-mediated synaptic vesicle fusion (PubMed:19386896). Required for the maintenance of synaptic ultrastructure in the adult retina (PubMed:19386896). Positively regulates synaptic transmission through synaptic vesicle availability and exocytosis of neurotransmitters at photoreceptor ribbon synapses in the retina (PubMed:15911881, PubMed:19386896, PubMed:27335398). Suppresses tonic photoreceptor activity and baseline 'noise' by suppression of Ca(2+) vesicle tonic release and the facilitation of evoked synchronous and asynchronous Ca(2+) vesicle release (PubMed:22694764, PubMed:27335398).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).