

Product datasheet for **TR507300**

Adamts19 Mouse shRNA Plasmid (Locus ID 240322)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Adamts19 Mouse shRNA Plasmid (Locus ID 240322)
Locus ID:	240322
Synonyms:	4831407I23Rik; AU015154; D230034E10Rik
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Adamts19 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 240322). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_175506 , NM_175506.1 , NM_175506.2 , NM_175506.3 , NM_175506.4 , BC150736
UniProt ID:	P59509
Summary:	This gene encodes a member of "a disintegrin and metalloproteinase with thrombospondin motifs" (ADAMTS) family of multi-domain matrix-associated metalloendopeptidases that have diverse roles in tissue morphogenesis and pathophysiological remodeling, in inflammation and in vascular biology. This gene is predominantly expressed in the ovary with lower levels of expression observed in kidney, heart, skeletal muscle, lung and testis. The encoded preproprotein undergoes proteolytic processing to generate an active protease. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).