

## Product datasheet for **TR508116**

### Lpcat2 Mouse shRNA Plasmid (Locus ID 270084)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Lpcat2 Mouse shRNA Plasmid (Locus ID 270084)
Locus ID:	270084
Synonyms:	A330042H22; Ayt11; Ayt11a; lpafat1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Lpcat2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 270084). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_173014</a> , <a href="#">NM_001357374</a> , <a href="#">NM_001357375</a> , <a href="#">NM_173014.1</a> , <a href="#">BC027692</a> , <a href="#">BC144805</a> , <a href="#">BC144806</a> , <a href="#">BC150666</a> , <a href="#">BC160287</a> , <a href="#">NM_173014.2</a>
UniProt ID:	<a href="#">Q8BYI6</a>
Summary:	Possesses both acyltransferase and acetyltransferase activities (PubMed:17182612, PubMed:18156367). Activity is calcium-dependent (PubMed:17182612). Involved in platelet-activating factor (PAF) biosynthesis by catalyzing the conversion of the PAF precursor, 1-O-alkyl-sn-glycero-3-phosphocholine (lyso-PAF) into 1-O-alkyl-2-acetyl-sn-glycero-3-phosphocholine (PAF) (PubMed:17182612). Also converts lyso-PAF to 1-O-alkyl-2-acyl-sn-glycero-3-phosphocholine (PC), a major component of cell membranes and a PAF precursor (PubMed:17182612, PubMed:18156367). Under resting conditions, acyltransferase activity is preferred (PubMed:17182612). Upon acute inflammatory stimulus, acetyltransferase activity is enhanced and PAF synthesis increases (PubMed:17182612). Also catalyzes the conversion of 1-acyl-sn-glycero-3-phosphocholine to 1,2-diacyl-sn-glycero-3-phosphocholine. Involved in the regulation of lipid droplet number and size (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).