

Product datasheet for **TR508547**

Col27a1 Mouse shRNA Plasmid (Locus ID 373864)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Col27a1 Mouse shRNA Plasmid (Locus ID 373864)
Locus ID:	373864
Synonyms:	5730512J02Rik; AI449266
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Col27a1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 373864). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_025685 , NM_025685.2 , BC038041 , BC039222 , BC049173
UniProt ID:	Q5QNO9
Summary:	This gene encodes the alpha-1 subunit of type XXVII collagen, one of the low abundance fibril-forming collagens found in cartilage. The encoded protein forms a homotrimeric triple helical procollagen that undergoes proteolytic processing during fibril formation. Transgenic mice lacking a portion of the collagenous domain in the encoded protein exhibit skeletal abnormalities, chondrodysplasia and die at birth because of a lung defect. [provided by RefSeq, Dec 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).