

Product datasheet for **TR509197**

Retsat Mouse shRNA Plasmid (Locus ID 67442)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Retsat Mouse shRNA Plasmid (Locus ID 67442)
Locus ID:	67442
Synonyms:	0610039N19Rik; C80029; MMT-7; Ppsig
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Retsat - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 67442). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC117751 , NM_026159 , NM_026159.1 , NM_026159.2 , NM_026159.3 , NM_026159.4 , BC011203 , BC117750 , NM_026159.5
UniProt ID:	Q64FW2
Summary:	Catalyzes the saturation of all-trans-retinol to all-trans-13,14-dihydroretinol (PubMed:15358783, PubMed:17253779, PubMed:19139408). Does not exhibit any activity toward all-trans-retinoic acid, nor 9-cis, 11-cis or 13-cis-retinol isomers (PubMed:15358783). May play a role in the metabolism of vitamin A (PubMed:15358783, PubMed:17253779). Independently of retinol conversion, may regulate liver metabolism upstream of MLXIPL/ChREBP (PubMed:28855500). Required for adipocyte differentiation in a 3T3-L1 cell culture model (PubMed:19139408). This effect seems not to mimic the in vivo situation in which animals show increased adiposity in the absence of RETSAT (PubMed:19940255). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).