

## Product datasheet for **TR509322**

### Myc Mouse shRNA Plasmid (Locus ID 17869)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Myc Mouse shRNA Plasmid (Locus ID 17869)
Locus ID:	17869
Synonyms:	AU016757; bHLHe3; bHLHe39; Myc2; N; Niard; Nird
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Myc - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 17869). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001177352</a> , <a href="#">NM_001177353</a> , <a href="#">NM_001177354</a> , <a href="#">NM_010849</a> , <a href="#">NM_010849.1</a> , <a href="#">NM_010849.2</a> , <a href="#">NM_010849.3</a> , <a href="#">NM_010849.4</a> , <a href="#">NM_001177352.1</a> , <a href="#">NM_001177354.1</a> , <a href="#">BC138931</a> , <a href="#">BC006728</a> , <a href="#">BC138932</a>
UniProt ID:	<a href="#">P01108</a>
Summary:	The protein encoded by this gene is a multifunctional, nuclear phosphoprotein that plays a role in cell cycle progression, apoptosis and cellular transformation. It functions as a transcription factor that regulates transcription of specific target genes. Mutations, overexpression, rearrangement and translocation of this gene have been associated with a variety of hematopoietic tumors, leukemias and lymphomas, including Burkitt lymphoma, in human. There is evidence to show that alternative translation initiations from an upstream, in-frame non-AUG (CUG) and a downstream AUG start site result in the production of two isoforms with distinct N-termini, in human and mouse. Under conditions of stress, such as high cell densities and methionine deprivation, there is a specific and dramatic increase in the synthesis of the non-AUG initiated protein, suggesting its importance in times of adversity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).