

Product datasheet for **TR510109**

Trp53bp2 Mouse shRNA Plasmid (Locus ID 209456)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Trp53bp2 Mouse shRNA Plasmid (Locus ID 209456)
Locus ID:	209456
Synonyms:	53BP2; AI746547; ASPP2; PPP1R13A; Tp53bp2; X98550
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Trp53bp2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 209456). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC042874 , NM_173378 , NM_173378.1 , BC025187 , BC030894
UniProt ID:	Q8CG79
Summary:	Regulator that plays a central role in regulation of apoptosis and cell growth via its interactions. Regulates p53/TP53 by enhancing the DNA binding and transactivation function of p53/TP53 on the promoters of proapoptotic genes in vivo. Inhibits the ability of APPBP1 to conjugate NEDD8 to CUL1, and thereby decreases APPBP1 ability to induce apoptosis. Impedes cell cycle progression at G2/M. Its apoptosis-stimulating activity is inhibited by its interaction with DDX42 (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).