

## Product datasheet for **TR510188**

### Aire Mouse shRNA Plasmid (Locus ID 11634)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Aire Mouse shRNA Plasmid (Locus ID 11634)
Locus ID:	11634
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Aire - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 11634). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC103510</a> , <a href="#">BC103511</a> , <a href="#">BC103518</a> , <a href="#">NM_001271549</a> , <a href="#">NM_001271550</a> , <a href="#">NM_001271551</a> , <a href="#">NM_001271552</a> , <a href="#">NM_001271553</a> , <a href="#">NM_001271554</a> , <a href="#">NM_001271555</a> , <a href="#">NM_001271556</a> , <a href="#">NM_001271557</a> , <a href="#">NM_001271558</a> , <a href="#">NM_001271559</a> , <a href="#">NM_009646</a> , <a href="#">NR_073358</a> , <a href="#">NM_009646.1</a> , <a href="#">NM_009646.2</a> , <a href="#">NM_001271559.1</a> , <a href="#">NM_001271558.1</a> , <a href="#">NM_001271557.1</a> , <a href="#">NM_001271556.1</a> , <a href="#">NM_001271555.1</a> , <a href="#">NM_001271554.1</a> , <a href="#">NM_001271553.1</a> , <a href="#">NM_001271552.1</a> , <a href="#">NM_001271551.1</a> , <a href="#">NM_001271550.1</a> , <a href="#">NM_001271549.1</a>
UniProt ID:	<a href="#">Q9Z0E3</a>



[View online »](#)

**Summary:** Transcription factor playing an essential role to promote self-tolerance in the thymus by regulating the expression of a wide array of self-antigens that have the commonality of being tissue-restricted in their expression pattern in the periphery, called tissue restricted antigens (TRA) (Probable). Binds to G-doublets in an A/T-rich environment; the preferred motif is a tandem repeat of 5'-ATTGGTTA-3' combined with a 5'-TTATTA-3' box. Binds to nucleosomes (By similarity). Binds to chromatin and interacts selectively with histone H3 that is not methylated at 'Lys-4', not phosphorylated at 'Thr-3' and not methylated at 'Arg-2'. Functions as a sensor of histone H3 modifications that are important for the epigenetic regulation of gene expression. Mainly expressed by medullary thymic epithelial cells (mTECs), induces the expression of thousands of tissue-restricted proteins, which are presented on major histocompatibility complex class I (MHC-I) and MHC-II molecules to developing T-cells percolating through the thymic medulla (By similarity). Also induces self-tolerance through other mechanisms such as the regulation of the mTEC differentiation program (PubMed:19015306). Controls the medullary accumulation of thymic dendritic cells and the development of regulatory T-cell through the regulation of XCL1 expression (PubMed:21300913). Regulates the production of CCR4 and CCR7 ligands in medullary thymic epithelial cells and alters the coordinated maturation and migration of thymocytes (PubMed:19923453). In thymic B-cells, allows the presentation of licensing-dependent endogenous self-antigen for negative selection (PubMed:26070482). In secondary lymphoid organs, induces functional inactivation of CD4(+) T-cells. Expressed by a distinct bone marrow-derived population, induces self-tolerance through a mechanism that does not require regulatory T-cells and is resistant to innate inflammatory stimuli (PubMed:23993652). [UniProtKB/Swiss-Prot Function]

**shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).