

## Product datasheet for **TR510216**

### **Rrp1b Mouse shRNA Plasmid (Locus ID 72462)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	Rrp1b Mouse shRNA Plasmid (Locus ID 72462)
Locus ID:	72462
Synonyms:	2600005C20Rik; D030064A17; Kiaa0179; mKIAA0179
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Rrp1b - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 72462). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC016569</a> , <a href="#">NM_001163734</a> , <a href="#">NM_028244</a> , <a href="#">NM_001163734.1</a> , <a href="#">NM_028244.1</a> , <a href="#">NM_028244.2</a>
UniProt ID:	<a href="#">Q91YK2</a>
Summary:	Positively regulates DNA damage-induced apoptosis by acting as a transcriptional coactivator of proapoptotic target genes of the transcriptional activator E2F1 (By similarity). Likely to play a role in ribosome biogenesis by targeting serine/threonine protein phosphatase PP1 to the nucleolus (By similarity). Involved in regulation of mRNA splicing (PubMed:23604122). Inhibits SIPA1 GTPase activity (PubMed:18081427). Involved in regulating expression of extracellular matrix genes (PubMed:18081427). Associates with chromatin and may play a role in modulating chromatin structure (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).