

Product datasheet for **TR510296**

Gimap4 Mouse shRNA Plasmid (Locus ID 107526)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Gimap4 Mouse shRNA Plasmid (Locus ID 107526)
Locus ID:	107526
Synonyms:	AU019574; E430007K16Rik; I; Ian-1; Ian1; Imap4
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Gimap4 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 107526). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC005577 , NM_001243199 , NM_001243200 , NM_001243201 , NM_174990 , NM_175048 , NM_174990.1 , NM_174990.2 , NM_174990.3 , NM_174990.4 , NM_175048.1 , NM_175048.2 , NM_175048.3 , NM_001243201.1 , NM_001243200.1 , NM_001243199.1 , BC098226
UniProt ID:	Q99JY3
Summary:	This gene encodes a protein belonging to the GTP-binding superfamily and to the immun-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. This gene exists within a cluster of other related genes located on mouse chromosome 6. This family member encodes a lymphoid signaling protein that functions to accelerate programmed T-cell death, which appears to correlate with the phosphorylation status of the protein. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).