

Product datasheet for **TR510319**

Nlrp1b Mouse shRNA Plasmid (Locus ID 637515)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Nlrp1b Mouse shRNA Plasmid (Locus ID 637515)
Locus ID:	637515
Synonyms:	ENSMUSG00000070390; Nalp1b; OTTMUSG00000006087; OTTMUSG00000006089
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Nlrp1b - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 637515). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001040696 , NM_001162414 , NM_001040696.1 , NM_001162414.1 , BC141354 , NM_001039680
UniProt ID:	A1Z198
Summary:	As the sensor component of the NLRP1 inflammasome, plays a crucial role in innate immunity and inflammation. In response to pathogens and other damage-associated signals, initiates the formation of the inflammasome polymeric complex, made of Nlrp1b, CASP1, and possibly PYCARD. Recruitment of proCASP1 to the inflammasome promotes its activation and CASP1-catalyzed IL1B and IL18 maturation and secretion in the extracellular milieu. Activation of NLRP1 inflammasome is also required for HMGB1 secretion. The active cytokines and HMGB1 stimulate inflammatory responses. Inflammasomes can also induce pyroptosis, an inflammatory form of programmed cell death (By similarity). May be activated by muramyl dipeptide (MDP), a fragment of bacterial peptidoglycan, in a NOD2-dependent manner (PubMed:18511561). Might be activated by Toxoplasma gondii, although at a lower extent than allele 1 (PubMed:24218483). Contrary to Nlrp1b allele 1, allele 2 is not activated by Bacillus anthracis lethal toxin (PubMed:16429160, PubMed:21170303, PubMed:24492532). [UniProtKB/Swiss-Prot Function]


[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>