

Product datasheet for **TR510416**

Ciita Mouse shRNA Plasmid (Locus ID 12265)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Ciita Mouse shRNA Plasmid (Locus ID 12265)
Locus ID:	12265
Synonyms:	C2t; C2ta; EG669998; Gm9475; Mhc2ta
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Ciita - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 12265). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001243760 , NM_001243761 , NM_001302618 , NM_001302619 , NM_007575 , NM_007575.1 , NM_007575.2 , NM_007575.3 , NM_007575.4 , NM_001302619.1 , NM_001243761.1 , NM_001243761.2 , NM_001243760.1 , NM_001243760.2 , NM_001302618.1 , BC107330 , BC107331 , BC140358 , BC148699
UniProt ID:	P79621
Summary:	This gene encodes a member of the NOD-like receptor protein family. This protein acts as a transcriptional coactivator and component of the enhanceosome complex to stimulate transcription of MHC class II genes in the adaptive immune response. This protein may also regulate the transcription of MHC class I genes. Mutations in the human gene have been linked to a rare immunodeficiency, bare lymphocyte syndrome, and homozygous knockout mice exhibit many features of this disease. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Oct 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).