

Product datasheet for **TR510502**

Runx2 Mouse shRNA Plasmid (Locus ID 12393)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Runx2 Mouse shRNA Plasmid (Locus ID 12393)
Locus ID:	12393
Synonyms:	AM; AML3; Cbf; Cbfa; Cbfa-1; Cbfa1; LS3; Os; Osf2; PEB; Pebp2a1; PEBP2aA; Pebpa2a
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Runx2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 12393). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001145920 , NM_001146038 , NM_001271627 , NM_001271630 , NM_001271631 , NM_009820 , NR_073392 , NR_073425 , NM_001146038.1 , NM_001146038.2 , NM_009820.1 , NM_009820.2 , NM_009820.3 , NM_009820.4 , NM_009820.5 , NM_001145920.1 , NM_001145920.2 , NM_001271631.1 , NM_001271630.1 , NM_001271627.1 , BC172607 , BC172715 , NM_001271633
UniProt ID:	Q08775
Summary:	This gene encodes a member of the runt domain-containing family of transcription factors. This protein is essential for osteoblastic differentiation and skeletal morphogenesis and acts as a scaffold for nucleic acids and regulatory factors involved in skeletal gene expression. The protein can bind DNA both as a monomer or, with more affinity, as a subunit of a heterodimeric complex. Transcript variants that encode different protein isoforms result from the use of alternate promoters as well as alternate splicing. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).