

Product datasheet for **TR510515**

Ercc8 Mouse shRNA Plasmid (Locus ID 71991)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Ercc8 Mouse shRNA Plasmid (Locus ID 71991)
Locus ID:	71991
Synonyms:	2410022P04Rik; 2810431L23Rik; 4631412O06Rik; B130065P18Rik; Ckn1; Csa
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Ercc8 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 71991). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC037200 , NM_028042 , NM_028042.1 , NM_028042.2 , NM_028042.3 , BC056463 , BC066170 , NM_001362403
UniProt ID:	Q8CFD5
Summary:	Involved in transcription-coupled nucleotide excision repair. It is required for the recruitment of XAB2, HMG1 and TCEA1/TFIIIS to a transcription-coupled repair complex which removes RNA polymerase II-blocking lesions from the transcribed strand of active genes. It is the substrate-recognition component of the CSA complex (DCX(ERCC8) complex) which promotes the ubiquitination and subsequent proteasomal degradation of ERCC6 in a UV-dependent manner; ERCC6 degradation is essential for the recovery of RNA synthesis after transcription-coupled repair.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).