

Product datasheet for **TR510822**

Kars Mouse shRNA Plasmid (Locus ID 85305)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Kars Mouse shRNA Plasmid (Locus ID 85305)
Locus ID:	85305
Synonyms:	AA589550; AL024334; AL033315; AL033367; D8Ertd698e; D8Wsu108e; LysRS; mKIAA0070
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Kars - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 85305). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC027356 , BC035324 , BC036289 , NM_001130868 , NM_001286384 , NM_053092 , NM_001130868.1 , NM_001130868.2 , NM_053092.1 , NM_053092.2 , NM_053092.3 , NM_001286384.1 , BC027356.1 , NM_001363429 , NM_001363430
UniProt ID:	Q99MN1
Summary:	Catalyzes the specific attachment of an amino acid to its cognate tRNA in a 2 step reaction: the amino acid (AA) is first activated by ATP to form AA-AMP and then transferred to the acceptor end of the tRNA. When secreted, acts as a signaling molecule that induces immune response through the activation of monocyte/macrophages. Catalyzes the synthesis of the signaling molecule diadenosine tetraphosphate (Ap4A), and thereby mediates disruption of the complex between HINT1 and MITF and the concomitant activation of MITF transcriptional activity.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).