

Product datasheet for **TR511027**

Foxn1 Mouse shRNA Plasmid (Locus ID 15218)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Foxn1 Mouse shRNA Plasmid (Locus ID 15218)
Locus ID:	15218
Synonyms:	D11Bhm185e; Fkh19; HFH-11; Hfh11; nu; nude; wh; Whn
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Foxn1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 15218). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC108980 , BC108981 , NM_001277290 , NM_008238 , NM_008238.1 , NM_008238.2 , NM_001277290.1
UniProt ID:	Q61575
Summary:	The protein encoded by this gene is part of the forkhead family or "winged-helix" transcription factors that are important in developmental processes, immune system regulation, metabolism, cancer and aging. This gene family has over 100 members, subdivided into classes (A-Q) based on phylogeny. The encoded protein is proposed to regulate development of the thymus and differentiation of keratinocytes. Mutations in this gene cause severe primary T-cell immunodeficiency and congenital alopecia. In mouse mutations of this gene underlie the phenotype of the nude mouse, which has been widely used as a model system in oncology, immunology, dermatology, and transplantation studies. In humans mutations in this gene have been correlated with T-cell immunodeficiency, the skin disorder congenital alopecia, and nail dystrophy. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Apr 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).