

Product datasheet for **TR511536**

Ptges2 Mouse shRNA Plasmid (Locus ID 96979)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Ptges2 Mouse shRNA Plasmid (Locus ID 96979)
Locus ID:	96979
Synonyms:	0610038H10Rik; C79137; Gbf1; Mptes2; Pges2
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Ptges2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 96979). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC004846 , NM_133783 , NM_133783.1 , NM_133783.2
UniProt ID:	Q8BWM0
Summary:	The protein encoded by this gene is a Golgi membrane-associated prostaglandin E synthase candidate, which is capable of catalyzing the conversion of prostaglandin H2 to prostaglandin E2 in vitro. However, a study using mice deficient of this gene suggests that this enzyme does not contribute to prostaglandin E2 biosynthesis in vivo. This protein is synthesized as a Golgi membrane-bound protein, but its N-terminal hydrophobic region is cleaved off during protein maturation to produce the predominant soluble truncated form that still retains the enzyme activity. This soluble protein also has been shown to activate the transcription regulated by a gamma-interferon-activated transcription element (GATE), possibly via an interaction with CAAAT/enhancer-binding protein-beta. [provided by RefSeq, Oct 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).