

Product datasheet for **TR512053**

Zfp451 Mouse shRNA Plasmid (Locus ID 98403)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Zfp451 Mouse shRNA Plasmid (Locus ID 98403)
Locus ID:	98403
Synonyms:	4930515K21Rik; 4933435G09Rik; AI596398; COASTER; Kiaa0576-hp; mKIAA1702; Znf451
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Zfp451 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 98403). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC019546 , BC062154 , NM_001290699 , NM_001290700 , NM_133817 , NM_001359274 , NM_001359275 , NM_133817.1 , NM_133817.2 , NM_133817.3 , NM_001290700.1 , NM_001290699.1 , BC062154.1 , BC022767 , BC024435 , BC026587 , BC043685
UniProt ID:	Q8C0P7
Summary:	E3 SUMO-protein ligase; has a preference for SUMO2 and SUMO3 and facilitates UBE2I/UBC9-mediated sumoylation of target proteins. Plays a role in protein SUMO2 modification in response to stress caused by DNA damage and by proteasome inhibitors (in vitro). Required for MCM4 sumoylation. Has no activity with SUMO1 (PubMed:26524493). Preferentially transfers an additional SUMO2 chain onto the SUMO2 consensus site 'Lys-11'. Negatively regulates transcriptional activation mediated by the SMAD4 complex in response to TGF-beta signaling. Inhibits EP300-mediated acetylation of histone H3 at 'Lys-9'. Plays a role in regulating the transcription of AR targets (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).