

Product datasheet for **TR512285**

Pard6a Mouse shRNA Plasmid (Locus ID 56513)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Pard6a Mouse shRNA Plasmid (Locus ID 56513)
Locus ID:	56513
Synonyms:	0710008C04Rik; 2610010A15Rik; Par-6; PAR-6A; Par6; PAR6alpha; Par6c; TAX40; Tip-40
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Pard6a - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 56513). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC049593 , BC061114 , NM_001047435 , NM_001047436 , NM_001286344 , NM_001286345 , NM_019695 , NM_001047435.1 , NM_001047435.2 , NM_019695.1 , NM_019695.2 , NM_019695.3 , NM_001047436.1 , NM_001047436.2 , NM_001286344.1 , NM_001286345.1
UniProt ID:	Q9Z101
Summary:	Adapter protein involved in asymmetrical cell division and cell polarization processes. Probably involved in the formation of epithelial tight junctions. Association with PARD3 may prevent the interaction of PARD3 with F11R/JAM1, thereby preventing tight junction assembly. The PARD6-PARD3 complex links GTP-bound Rho small GTPases to atypical protein kinase C proteins (PubMed:15761148). Regulates centrosome organization and function. Essential for the centrosomal recruitment of key proteins that control centrosomal microtubule organization (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).