

Product datasheet for **TR512476**

Cd200r3 Mouse shRNA Plasmid (Locus ID 74603)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Cd200r3 Mouse shRNA Plasmid (Locus ID 74603)
Locus ID:	74603
Synonyms:	4733401118Rik; 4833409J19Rik; AI505817; BB106877; mCD200RLb
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Cd200r3 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 74603). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC106838 , NM_001128132 , NM_001128133 , NM_027578 , NM_029018 , NM_027578.1 , NM_001128132.1 , NM_029018.1 , NM_029018.2 , NM_029018.3 , NM_029018.4 , NM_001128133.1
UniProt ID:	Q5UKY4
Summary:	According to PubMed:15187158 isoform 4 is a receptor for the CD200 cell surface glycoprotein. According to PubMed:16081818 isoform 4 is not a receptor for the CD200/OX2 cell surface glycoprotein. Isoform 1, isoform 2 and isoform 3 are involved in the recruitment or surface expression of the TYROBP receptor. Isoform 6, isoform 7 and isoform 8 are not involved in the recruitment or surface expression of the TYROBP receptor.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).