

## Product datasheet for **TR512513**

### Plk3 Mouse shRNA Plasmid (Locus ID 12795)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Plk3 Mouse shRNA Plasmid (Locus ID 12795)
Locus ID:	12795
Synonyms:	Cn; Cnk; Fnk; PLK-3; PRK
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Plk3 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 12795). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC031180</a> , <a href="#">BC063051</a> , <a href="#">NM_013807</a> , <a href="#">NM_013807.1</a> , <a href="#">NM_013807.2</a> , <a href="#">NM_013807.3</a>
Summary:	This gene encodes a member of the highly conserved polo-like kinase family of serine/threonine kinases. Members of this family are characterized by an amino-terminal catalytic domain and a carboxy-terminal bipartite polo box domain that functions as a substrate-binding motif and a cellular localization signal. Polo-like kinases have primarily been implicated in cell cycle regulation. In mouse, this protein that has been reported to localize to the nucleolus during interphase but is undetectable during mitosis, following nucleolus dissociation during prophase. The protein relocates to the nucleolus just prior to cytokinesis and peak levels are detected during G1 of interphase. This gene has been implicated in regulation of entry into S phase, with RNAi-induced depletion resulting in failure to re-enter the cell cycle. Mice deficient for this gene exhibit increased weight and tumor development at advanced age. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).