

Product datasheet for **TR512643**

Smchd1 Mouse shRNA Plasmid (Locus ID 74355)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Smchd1 Mouse shRNA Plasmid (Locus ID 74355)
Locus ID:	74355
Synonyms:	4931400A14Rik; AW554188; mKIAA0650; MommeD1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Smchd1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 74355). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC062946 , NM_028887 , NM_028887.1 , NM_028887.2 , NM_028887.3 , BC039785 , BC040778 , BC044905 , BC058205 , BC058618 , BC076627
UniProt ID:	Q6P5D8



[View online »](#)

Summary: Non-canonical member of the structural maintenance of chromosomes (SMC) protein family that plays a key role in epigenetic silencing by regulating chromatin architecture (PubMed:26091879, PubMed:29887375). Promotes heterochromatin formation in both autosomes and chromosome X, probably by mediating the merge of chromatin compartments (PubMed:23754746, PubMed:23819640, PubMed:26391951, PubMed:28587678, PubMed:29887375). Plays a key role in chromosome X inactivation in females by promoting the spreading of heterochromatin (PubMed:18425126, PubMed:22841499, PubMed:26391951, PubMed:29887375). Recruited to inactivated chromosome X by Xist RNA and acts by mediating the merge of chromatin compartments: promotes random chromatin interactions that span the boundaries of existing structures, leading to create a compartment-less architecture typical of inactivated chromosome X (PubMed:29887375). Required to facilitate Xist RNA spreading (PubMed:29887375). Also required for silencing of a subset of clustered autosomal loci in somatic cells, such as the DUX4 locus (PubMed:23754746, PubMed:23819640, PubMed:28587678). Has ATPase activity; may participate in structural manipulation of chromatin in an ATP-dependent manner as part of its role in gene expression regulation (PubMed:26391951, PubMed:27059856). Also plays a role in DNA repair: localizes to sites of DNA double-strand breaks in response to DNA damage to promote the repair of DNA double-strand breaks (By similarity). Acts by promoting non-homologous end joining (NHEJ) and inhibiting homologous recombination (HR) repair (By similarity). Required during preimplantation development, probably acts by regulating chromatin architecture (PubMed:29900695).[UniProtKB/Swiss-Prot Function]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).