

## Product datasheet for **TR512872**

### Lmnb2 Mouse shRNA Plasmid (Locus ID 16907)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Lmnb2 Mouse shRNA Plasmid (Locus ID 16907)
Locus ID:	16907
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Lmnb2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 16907). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC042430</a> , <a href="#">BC051985</a> , <a href="#">NM_010722</a> , <a href="#">NM_010722.2</a> , <a href="#">NM_010722.3</a> , <a href="#">NM_010722.4</a> , <a href="#">NM_010722.5</a>
UniProt ID:	<a href="#">P21619</a>
Summary:	This gene encodes a protein component of the nuclear lamina, which provides a structural framework for the nuclear envelope. Defects in this gene were found to cause abnormalities in the shape of neurons. This locus represents one of two B-type lamin genes that may be partially, but not entirely, functionally redundant in neuronal development. Loss of both B-type lamin genes in keratinocytes results in ichthyosis and a skin barrier defect leading to dehydration. Alternative transcriptional initiation and splicing results in multiple transcript variants and protein isoforms, including an isoform with a shorter N-terminal rod domain that may function in nuclear envelope remodeling during spermatogenesis. A related pseudogene is found on chromosome 5. [provided by RefSeq, Sep 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).