

Product datasheet for **TR513467**

Zfyve27 Mouse shRNA Plasmid (Locus ID 319740)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Zfyve27 Mouse shRNA Plasmid (Locus ID 319740)
Locus ID:	319740
Synonyms:	2210011N02Rik; 9530077C24Rik; AI426636; AI593546; AI835681
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Zfyve27 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 319740). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC042595 , NM_001164531 , NM_177319 , NM_177319.1 , NM_177319.2 , NM_177319.3 , NM_001164531.1 , BC024801
UniProt ID:	Q3TXX3
Summary:	Key regulator of RAB11-dependent vesicular trafficking during neurite extension through polarized membrane transport (By similarity). Promotes axonal elongation and contributes to the establishment of neuronal cell polarity (PubMed:24251978). Involved in nerve growth factor-induced neurite formation in VAPA-dependent manner. Contributes to both the formation and stabilization of the tubular ER network. Involved in ER morphogenesis by regulating the sheet-to-tubule balance and possibly the density of tubule interconnections (By similarity). Acts as an adapter protein that facilitates the interaction of KIF5A with VAPA, VAPB, SURF4, RAB11A, RAB11B and RTN3 and the ZFYVE27-KIF5A complex contributes to the transport of these proteins in neurons. Can induce formation of neurite-like membrane protrusions in non-neuronal cells in a KIF5A/B-dependent manner (PubMed:21976701). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).