

## Product datasheet for **TR513587**

### Mbnl2 Mouse shRNA Plasmid (Locus ID 105559)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Mbnl2 Mouse shRNA Plasmid (Locus ID 105559)
Locus ID:	105559
Synonyms:	1110002M11Rik; AI047808; AI837313; AI849185; AL118326; MBLL; mKIAA4072; R75232
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Mbnl2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 105559). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC075665</a> , <a href="#">NM_175341</a> , <a href="#">NM_207515</a> , <a href="#">NM_001360376</a> , <a href="#">NM_001360377</a> , <a href="#">NM_001360378</a> , <a href="#">NM_001360379</a> , <a href="#">NM_001360380</a> , <a href="#">NM_175341.1</a> , <a href="#">NM_175341.2</a> , <a href="#">NM_175341.3</a> , <a href="#">NM_175341.4</a> , <a href="#">NM_207515.1</a> , <a href="#">NM_207515.2</a> , <a href="#">BC020406</a> , <a href="#">BC023895</a> , <a href="#">BC026711</a> , <a href="#">BC079898</a> , <a href="#">BM932374</a>
UniProt ID:	<a href="#">Q8C181</a>
Summary:	Mediates pre-mRNA alternative splicing regulation. Acts either as activator or repressor of splicing on specific pre-mRNA targets. Inhibits cardiac troponin-T (TNNT2) pre-mRNA exon inclusion but induces insulin receptor (IR) pre-mRNA exon inclusion in muscle. Antagonizes the alternative splicing activity pattern of CELF proteins. RNA-binding protein that binds to 5'ACACCC-3' core sequence, termed zipcode, within the 3' UTR of ITGA3. Binds to CUG triplet repeat expansion in myotonic dystrophy muscle cells by sequestering the target RNAs. Seems to regulate expression and localization of ITGA3 by transporting it from the nucleus to cytoplasm at adhesion plaques. May play a role in myotonic dystrophy pathophysiology (DM) (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).