

Product datasheet for **TR513709**

Zfp598 Mouse shRNA Plasmid (Locus ID 213753)

Product data:

| | |
|---------------------------|--|
| Product Type: | shRNA Plasmids |
| Product Name: | Zfp598 Mouse shRNA Plasmid (Locus ID 213753) |
| Locus ID: | 213753 |
| Synonyms: | BC023040; Ntrap; Znf598 |
| Vector: | pRS (TR20003) |
| E. coli Selection: | Ampicillin |
| Mammalian Cell Selection: | Puromycin |
| Format: | Retroviral plasmids |
| Components: | Zfp598 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 213753). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free. |
| RefSeq: | BC050859 , NM_183149 , NM_001348231 , NM_183149.1 , NM_183149.2 , BC023040 , BC024690 , BC038138 |
| UniProt ID: | Q80YR4 |
| Summary: | E3 ubiquitin-protein ligase that plays a key role in the ribosome quality control (RQC), a pathway that takes place when a ribosome has stalled during translation. Required for ribosomes to terminally stall during translation of poly(A) sequences by mediating monoubiquitination of 40S ribosomal protein RPS10/eS10, RPS20/uS10 and RPS3/uS3. Stalling precludes synthesis of a long poly-lysine tail and initiates the RQC pathway to degrade the potentially detrimental aberrant nascent polypeptide. Also acts as a component of the 4EHP-GYF2 complex, a multiprotein complex that acts as a repressor of translation initiation.[UniProtKB/Swiss-Prot Function] |
| shRNA Design: | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service . |


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).