

## Product datasheet for **TR513768**

### Abi1 Mouse shRNA Plasmid (Locus ID 11308)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Abi1 Mouse shRNA Plasmid (Locus ID 11308)
Locus ID:	11308
Synonyms:	abi-1; E3B1; NAP1; Ssh3bp1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Abi1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 11308). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC004657</a> , <a href="#">BC079642</a> , <a href="#">NM_001077190</a> , <a href="#">NM_001077192</a> , <a href="#">NM_001077193</a> , <a href="#">NM_001331043</a> , <a href="#">NM_001331044</a> , <a href="#">NM_001331045</a> , <a href="#">NM_007380</a> , <a href="#">NM_145994</a> , <a href="#">NM_007380.1</a> , <a href="#">NM_007380.2</a> , <a href="#">NM_007380.3</a> , <a href="#">NM_007380.4</a> , <a href="#">NM_001077192.1</a> , <a href="#">NM_001077192.2</a> , <a href="#">NM_001077192.3</a> , <a href="#">NM_001077190.1</a> , <a href="#">NM_001077190.2</a> , <a href="#">NM_001077190.3</a> , <a href="#">NM_145994.1</a> , <a href="#">NM_145994.2</a> , <a href="#">NM_145994.3</a> , <a href="#">NM_001077193.1</a> , <a href="#">NM_001077193.2</a> , <a href="#">NM_001077193.3</a>
UniProt ID:	<a href="#">Q8CBW3</a>
Summary:	May act in negative regulation of cell growth and transformation by interacting with nonreceptor tyrosine kinases ABL1 and/or ABL2. In vitro, at least isoform 2 and isoform 4 suppress the transforming activity of Abelson murine leukemia virus (v-Abl) after overexpression in fibroblasts. May play a role in regulation EGF-induced Erk pathway activation. Involved in cytoskeletal reorganization and EGFR signaling. Together with EPS8 participates in transduction of signals from Ras to Rac. In vitro, a trimeric complex of ABI1, EPS8 and SOS1 exhibits Rac specific guanine nucleotide exchange factor (GEF) activity and ABI1 seems to act as an adapter in the complex. Regulates ABL1/c-Abl-mediated phosphorylation of ENAH. Recruits WASF1 to lamellipodia and there seems to regulate WASF1 protein level. In brain, seems to regulate the dendritic outgrowth and branching as well as to determine the shape and number of synaptic contacts of developing neurons. [UniProtKB/Swiss-Prot Function]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).