

## Product datasheet for **TR513847**

### Uimc1 Mouse shRNA Plasmid (Locus ID 20184)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Uimc1 Mouse shRNA Plasmid (Locus ID 20184)
Locus ID:	20184
Synonyms:	9430016E08Rik; D330018D10Rik; D630032M02Rik; RIP110; Rxrip110
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Uimc1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 20184). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC040808</a> , <a href="#">NM_001293660</a> , <a href="#">NM_011307</a> , <a href="#">NR_121592</a> , <a href="#">NM_001359757</a> , <a href="#">NM_001359758</a> , <a href="#">NM_011307.1</a> , <a href="#">NM_011307.2</a> , <a href="#">NM_001293660.1</a> , <a href="#">BC037092</a> , <a href="#">BC051659</a> , <a href="#">BC066220</a> , <a href="#">NM_011307.3</a>
UniProt ID:	<a href="#">Q5U5Q9</a>
Summary:	Ubiquitin-binding protein. Specifically recognizes and binds 'Lys-63'-linked ubiquitin (PubMed:19536136). Plays a central role in the BRCA1-A complex by specifically binding 'Lys-63'-linked ubiquitinated histones H2A and H2AX at DNA lesions sites, leading to target the BRCA1-BARD1 heterodimer to sites of DNA damage at double-strand breaks (DSBs). The BRCA1-A complex also possesses deubiquitinase activity that specifically removes 'Lys-63'-linked ubiquitin on histones H2A and H2AX. Also weakly binds monoubiquitin but with much less affinity than 'Lys-63'-linked ubiquitin. May interact with monoubiquitinated histones H2A and H2B; the relevance of such results is however unclear in vivo. Does not bind Lys-48'-linked ubiquitin. May indirectly act as a transcriptional repressor by inhibiting the interaction of NR6A1 with the corepressor NCOR1 (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).