

## Product datasheet for **TR513996**

### Cnot6l Mouse shRNA Plasmid (Locus ID 231464)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cnot6l Mouse shRNA Plasmid (Locus ID 231464)
Locus ID:	231464
Synonyms:	4932442K20Rik
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Cnot6l - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 231464). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC018506</a> , <a href="#">NM_001285511</a> , <a href="#">NM_001285514</a> , <a href="#">NM_144910</a> , <a href="#">NM_178854</a> , <a href="#">NM_144910.1</a> , <a href="#">NM_144910.2</a> , <a href="#">NM_178854.1</a> , <a href="#">NM_178854.2</a> , <a href="#">NM_178854.3</a> , <a href="#">NM_178854.4</a> , <a href="#">NM_001285514.1</a> , <a href="#">NM_001285511.1</a> , <a href="#">BC023923</a> , <a href="#">BC060137</a>
UniProt ID:	<a href="#">Q8VEG6</a>
Summary:	Poly(A) nuclease with 3'-5' RNase activity. Catalytic component of the CCR4-NOT complex which is one of the major cellular mRNA deadenylases and is linked to various cellular processes including bulk mRNA degradation, miRNA-mediated repression, translational repression during translational initiation and general transcription regulation. Additional complex functions may be a consequence of its influence on mRNA expression. Involved in mRNA decay mediated by the major-protein-coding determinant of instability (mCRD) of the FOS gene in the cytoplasm. Involved in deadenylation-dependent degradation of CDKN1B mRNA. Its mRNA deadenylase activity can be inhibited by TOB1. Mediates cell proliferation and cell survival and prevents cellular senescence (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).