

Product datasheet for **TR514868**

Abcg8 Mouse shRNA Plasmid (Locus ID 67470)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Abcg8 Mouse shRNA Plasmid (Locus ID 67470)
Locus ID:	67470
Synonyms:	1300003C16Rik; AI114946; sterolin-2
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Abcg8 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 67470). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001286005 , NM_001347418 , NM_026180 , NR_104382 , NM_026180.2 , NM_026180.3 , NM_001286005.1 , BC138484 , BC145293
UniProt ID:	Q9DBM0
Summary:	ABCG5 and ABCG8 form an obligate heterodimer that mediates Mg(2+)- and ATP-dependent sterol transport across the cell membrane (PubMed:16352607, PubMed:16867993, PubMed:18402465). Plays an essential role in the selective transport of the dietary cholesterol in and out of the enterocytes and in the selective sterol excretion by the liver into bile (PubMed:12444248, PubMed:14504269, PubMed:14657202, PubMed:25378657). Plays an important role in preventing the accumulation of dietary plant sterols in the body (PubMed:12444248, PubMed:14657202). Required for normal sterol homeostasis (PubMed:12444248, PubMed:14657202). The heterodimer with ABCG5 has ATPase activity (PubMed:16352607, PubMed:16867993).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).