

Product datasheet for **TR514955**

Kcnc2 Mouse shRNA Plasmid (Locus ID 268345)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Kcnc2 Mouse shRNA Plasmid (Locus ID 268345)
Locus ID:	268345
Synonyms:	AW047325; B230117I07; KShIIIA; Kv3.2
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Kcnc2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 268345). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<u>BC116289</u> , <u>BC116290</u> , <u>NM_001025581</u> , <u>NM_001359752</u> , <u>NM_001359753</u> , <u>NM_001025581.1</u>



[View online »](#)

Summary:

Voltage-gated potassium channel that mediates transmembrane potassium transport in excitable membranes, primarily in the brain. Contributes to the regulation of the fast action potential repolarization and in sustained high-frequency firing in neurons of the central nervous system (PubMed:10561420, PubMed:10414303, PubMed:11124984, PubMed:10903572, PubMed:11506885, PubMed:15317859, PubMed:15917463, PubMed:17761775, PubMed:21414897). Homotetramer channels mediate delayed-rectifier voltage-dependent potassium currents that activate rapidly at high-threshold voltages and inactivate slowly (PubMed:10414303). Forms tetrameric channels through which potassium ions pass in accordance with their electrochemical gradient. The channel alternates between opened and closed conformations in response to the voltage difference across the membrane (By similarity). Can form functional homotetrameric and heterotetrameric channels that contain variable proportions of KCNC1, and possibly other family members as well; channel properties depend on the type of alpha subunits that are part of the channel (PubMed:10531438, PubMed:12000114). Channel properties may be modulated by either the association with ancillary subunits, such as KCNE1, KCNE2 and KCNE3 or indirectly by nitric oxide (NO) through a cGMP- and PKG-mediated signaling cascade, slowing channel activation and deactivation of delayed rectifier potassium channels (By similarity). Contributes to fire sustained trains of very brief action potentials at high frequency in thalamocortical and suprachiasmatic nucleus (SCN) neurons, in hippocampal and neocortical interneurons and in retinal ganglion cells (PubMed:10561420, PubMed:10903572, PubMed:11506885, PubMed:17761775). Sustained maximal action potential firing frequency in inhibitory hippocampal interneurons is negatively modulated by histamine H2 receptor activation in a cAMP- and protein kinase (PKA) phosphorylation-dependent manner (PubMed:10903572). Plays a role in maintaining the fidelity of synaptic transmission in neocortical GABAergic interneurons by generating action potential (AP) repolarization at nerve terminals, thus reducing spike-evoked calcium influx and GABA neurotransmitter release (PubMed:15917463). Required for long-range synchronization of gamma oscillations over distance in the neocortex (PubMed:22539821). Contributes to the modulation of the circadian rhythm of spontaneous action potential firing in suprachiasmatic nucleus (SCN) neurons in a light-dependent manner (PubMed:21414897).[UniProtKB/Swiss-Prot Function]

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).