

## Product datasheet for **TR515366**

### Mib1 Mouse shRNA Plasmid (Locus ID 225164)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | Mib1 Mouse shRNA Plasmid (Locus ID 225164)   |
| Locus ID:                 | 225164   |
| Synonyms:                 | DIP-1; E430019M12Rik; Mib; mKIAA1323   |
| Vector:                   | pRS (TR20003)  |
| E. coli Selection:        | Ampicillin   |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Retroviral plasmids  |
| Components:               | Mib1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 225164). 5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.  |
| RefSeq:                   | <a href="#">BC131662</a> , <a href="#">NM_144860</a> , <a href="#">NM_144860.1</a> , <a href="#">NM_144860.2</a> , <a href="#">BC011287</a> , <a href="#">BC083072</a> , <a href="#">NM_001364997</a>  |
| UniProt ID:               | <a href="#">Q80SY4</a>   |
| Summary:                  | E3 ubiquitin-protein ligase that mediates ubiquitination of Delta receptors, which act as ligands of Notch proteins. Positively regulates the Delta-mediated Notch signaling by ubiquitinating the intracellular domain of Delta, leading to endocytosis of Delta receptors. Involved in ubiquitination of centriolar satellite CEP131, CEP290 and PCM1 proteins and hence inhibits primary cilium formation in proliferating cells. Mediates 'Lys-63'-linked polyubiquitination of TBK1, which probably participates in kinase activation (By similarity). Probably mediates ubiquitination and subsequent proteasomal degradation of DAPK1, thereby antagonizing anti-apoptotic effects of DAPK1 to promote TNF-induced apoptosis. [UniProtKB/Swiss-Prot Function] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).