

## Product datasheet for **TR516510**

### Jagn1 Mouse shRNA Plasmid (Locus ID 67767)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Jagn1 Mouse shRNA Plasmid (Locus ID 67767)
Locus ID:	67767
Synonyms:	5830427H10Rik
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Jagn1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 67767). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC011296</a> , <a href="#">BC025910</a> , <a href="#">BC037019</a> , <a href="#">NM_001205025</a> , <a href="#">NM_001243739</a> , <a href="#">NM_026365</a> , <a href="#">NR_037974</a> , <a href="#">NM_026365.1</a> , <a href="#">NM_026365.2</a> , <a href="#">NM_026365.3</a> , <a href="#">NM_001243739.1</a> , <a href="#">NM_001205025.1</a>
UniProt ID:	<a href="#">Q5XKN4</a>
Summary:	Endoplasmic reticulum transmembrane protein involved in vesicle-mediated transport, which is required for neutrophil function. Required for vesicle-mediated transport; it is however unclear whether it is involved in early secretory pathway or intracellular protein transport. Acts as a regulator of neutrophil function, probably via its role in vesicle-mediated transport: required for defense against fungal pathogens and for granulocyte colony-stimulating factor (GM-CSF) signaling pathway; possibly by regulating glycosylation and/or targeting of proteins contributing to the viability and migration of neutrophils.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).