

Product datasheet for **TR516905**

Otop1 Mouse shRNA Plasmid (Locus ID 21906)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Otop1 Mouse shRNA Plasmid (Locus ID 21906)
Locus ID:	21906
Synonyms:	A530025J20Rik; Otp1; tlt
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Otop1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 21906). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_172709 , NM_178139 , NM_172709.1 , NM_172709.2 , NM_172709.3 , BC146285 , BC148738
UniProt ID:	Q80VM9
Summary:	Proton-selective channel that specifically transports protons into cells (PubMed:29371428). Proton channel activity is only weakly-sensitive to voltage (PubMed:29371428). Proton-selective channel activity is probably required in cell types that use changes in intracellular pH for cell signaling or to regulate biochemical or developmental processes (PubMed:29371428). In the vestibular system of the inner ear, required for the formation and function of otoconia, which are calcium carbonate crystals that sense gravity and acceleration (PubMed:12651873). Probably acts by maintaining the pH appropriate for formation of otoconia (PubMed:29371428). Regulates purinergic control of intracellular calcium in vestibular supporting cells (PubMed:17606897, PubMed:20554841). May be involved in sour taste perception in sour taste cells by mediating entry of protons within the cytosol (PubMed:29371428). Also involved in energy metabolism, by reducing adipose tissue inflammation and protecting from obesity-induced metabolic dysfunction (PubMed:24379350).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).