

## Product datasheet for **TR516974**

### **Acsl1 Mouse shRNA Plasmid (Locus ID 14081)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	Acsl1 Mouse shRNA Plasmid (Locus ID 14081)
Locus ID:	14081
Synonyms:	Ac; Acas; Acas1; Acs; Fac; FacI2; FACS; LACS 1; LACS1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Acsl1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 14081). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC056644</a> , <a href="#">NM_001302163</a> , <a href="#">NM_007981</a> , <a href="#">NM_007981.1</a> , <a href="#">NM_007981.2</a> , <a href="#">NM_007981.3</a> , <a href="#">NM_007981.4</a> , <a href="#">NM_001302163.1</a> , <a href="#">BC006692</a>
UniProt ID:	<a href="#">P41216</a>
Summary:	The protein encoded by this gene belongs to a family of acyl coenzyme A synthetase proteins, which convert long chain fatty acids to acyl CoA products via an ATP-dependent pathway. This enzyme is enriched in heart, liver and adipose tissue, where it functions in lipid synthesis and mitochondrial and peroxisomal beta-oxidation. In addition, it is expressed in monocytes and macrophages where it appears to have a functionally distinct role in mediating inflammatory and innate immune responses. A pseudogene of this gene is found on chromosome 5. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).