

## Product datasheet for **TR517279**

### Tlr9 Mouse shRNA Plasmid (Locus ID 81897)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Tlr9 Mouse shRNA Plasmid (Locus ID 81897)
Locus ID:	81897
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Tlr9 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 81897). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_031178</a> , <a href="#">NM_031178.1</a> , <a href="#">NM_031178.2</a> , <a href="#">BC172644</a>
UniProt ID:	<a href="#">Q9EQU3</a>
Summary:	Key component of innate and adaptive immunity. TLRs (Toll-like receptors) control host immune response against pathogens through recognition of molecular patterns specific to microorganisms. TLR9 is a nucleotide-sensing TLR which is activated by unmethylated cytidine-phosphate-guanosine (CpG) dinucleotides. Acts via MYD88 and TRAF6, leading to NF-kappa-B activation, cytokine secretion and the inflammatory response (PubMed:18931679, PubMed:21402738, PubMed:14993594, PubMed:17474149, PubMed:25686612, PubMed:18820679). Plays a role in defense against systemic mouse cytomegalovirus infection (PubMed:14993594). Controls lymphocyte response to Helicobacter infection (PubMed:17474149). Upon CpG stimulation, induces B-cell proliferation, activation, survival and antibody production (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).