

## Product datasheet for **TR517675**

### **Rbm15 Mouse shRNA Plasmid (Locus ID 229700)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	Rbm15 Mouse shRNA Plasmid (Locus ID 229700)
Locus ID:	229700
Synonyms:	C230088J01Rik; mKIAA1438
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Rbm15 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 229700). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC120590</a> , <a href="#">NM_001045807</a> , <a href="#">NM_001045807.1</a> , <a href="#">BC028452</a> , <a href="#">BC051409</a> , <a href="#">BC057038</a> , <a href="#">BC080828</a> , <a href="#">BC099595</a> , <a href="#">BC137741</a>
UniProt ID:	<a href="#">Q0VBL3</a>



[View online »](#)

**Summary:**

RNA-binding protein that acts as a key regulator of N6-methyladenosine (m6A) methylation of RNAs, thereby regulating different processes, such as hematopoietic cell homeostasis, alternative splicing of mRNAs and X chromosome inactivation mediated by Xist RNA (PubMed:29535189). Associated component of the WMM complex, a complex that mediates N6-methyladenosine (m6A) methylation of RNAs, a modification that plays a role in the efficiency of mRNA splicing and RNA processing (PubMed:29535189). Plays a key role in m6A methylation, possibly by binding target RNAs and recruiting the WMM complex (PubMed:29535189). Involved in random X inactivation mediated by Xist RNA: acts by binding Xist RNA and recruiting the WMM complex, which mediates m6A methylation, leading to target YTHDC1 reader on Xist RNA and promoting transcription repression activity of Xist (By similarity). Required for the development of multiple tissues, such as the maintenance of the homeostasis of long-term hematopoietic stem cells and for megakaryocyte (MK) and B-cell differentiation (PubMed:17283045, PubMed:17376872, PubMed:18981216, PubMed:25468569). Regulates megakaryocyte differentiation by regulating alternative splicing of genes important for megakaryocyte differentiation; probably regulates alternative splicing via m6A regulation (By similarity). Required for placental vascular branching morphogenesis and embryonic development of the heart and spleen (PubMed:18981216). Acts as a regulator of thrombopoietin response in hematopoietic stem cells by regulating alternative splicing of MPL (PubMed:25468569). May also function as an mRNA export factor, stimulating export and expression of RTE-containing mRNAs which are present in many retrotransposons that require to be exported prior to splicing (By similarity). High affinity binding of pre-mRNA to RBM15 may allow targeting of the mRNP to the export helicase DBP5 in a manner that is independent of splicing-mediated NXF1 deposition, resulting in export prior to splicing (By similarity). May be implicated in HOX gene regulation (By similarity). [UniProtKB/Swiss-Prot Function]

**shRNA Design:**

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).