

Product datasheet for **TR518209**

Gm7609 Mouse shRNA Plasmid (Locus ID 665378)

Product data:

Product Type:	shRNA Plasmids
Locus ID:	665378
Synonyms:	EG665378
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Gm7609 – Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 665378). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001081746 , NM_001081746.1 , BC048373 , BC059276 , BC080715
Summary:	This gene is one of several duplicated genes in a repeat cluster in the homogeneously staining region (HSR) of chromosome 1. It is an Sp100-rs gene copy, which arose by duplication following the fusion of the Csprs (component of Sp100-rs) gene and the 5' part of the Sp100 (nuclear antigen Sp100) gene. The coding sequence of this gene is intact compared to that of the Csprs gene, and therefore NCBI represents this as a protein-coding gene. However, it is also possible to interpret this fusion gene as a pseudogene of the Sp100 gene. It is unknown if this gene produces a protein in vivo. [provided by RefSeq, Mar 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).