

Product datasheet for **TR518214**

Slurp2 Mouse shRNA Plasmid (Locus ID 69462)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Slurp2 Mouse shRNA Plasmid (Locus ID 69462)
Locus ID:	69462
Synonyms:	2300005B03Rik; SLURP-2
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Slurp2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 69462). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC115610 , BC115611 , NM_001081961 , NM_001081961.1
UniProt ID:	P0DP59
Summary:	Binds and may modulate the functional properties of nicotinic and muscarinic acetylcholine receptors. May regulate keratinocytes proliferation, differentiation and apoptosis. In vitro moderately inhibits ACh-evoked currents of alpha-3:beta-2-containing nAChRs, strongly these of alpha-4:beta-2-containing nAChRs, modulates alpha-7-containing nAChRs, and inhibits nicotine-induced signaling probably implicating alpha-3:beta-4-containing nAChRs. Proposed to act on alpha-3:beta-2 and alpha-7 nAChRs in an orthosteric, and on mAChRs, such as CHRM1 and CHRM3, in an allosteric manner.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).