

## Product datasheet for **TR518861**

### Sap18b Mouse shRNA Plasmid (Locus ID 100041953)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Sap18b Mouse shRNA Plasmid (Locus ID 100041953)
Locus ID:	100041953
Synonyms:	ENSMUSG00000061104; Gm10094
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Sap18b - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 100041953). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001142441</a> , <a href="#">NM_001142441.1</a>
Summary:	This intronless gene is highly similar to the multi-exon Sap18 gene on chromosome 14, whose product functions in transcriptional repression as a component of the Sin3 histone deacetylase complex. This gene may possibly be a Sap18 pseudogene, but it is represented as protein-coding because it appears to be transcribed and has an intact ORF that would result in a protein that is 100% identical to the Sap18 protein. [provided by RefSeq, Dec 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).