

Product datasheet for **TR519032**

Cdh10 Mouse shRNA Plasmid (Locus ID 320873)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Cdh10 Mouse shRNA Plasmid (Locus ID 320873)
Locus ID:	320873
Synonyms:	A830016G23Rik; C030003B10Rik; C030011H18Rik; T2-cadh
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Cdh10 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 320873). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC052056 , BC062962 , NM_009865 , NM_009865.1 , NM_009865.2 , NM_009865.3 , BC055808
UniProt ID:	P70408
Summary:	This gene encodes a member of the cadherin family of calcium-dependent glycoproteins that mediate cell adhesion and regulate many morphogenetic events during development. The encoded preproprotein is further processed to generate a mature protein. This gene is expressed in the blood-brain barrier and retinal endothelia suggesting a role in the development and maintenance of brain barrier. Alternative splicing results in multiple transcript variants. Multiple distinct genes of the cadherin family, including this gene, are found on chromosome 15. [provided by RefSeq, Oct 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).